




IKZF3 Antibody

Product Code	CSB-PA892336DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKT9
Immunogen	Recombinant Human Zinc finger protein Aiolos protein (280-450AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein Aiolos (Ikaros family zinc finger protein 3), IKZF3, ZNFN1A3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	IKZF3

Image



Western blot

All lanes: Zinc finger protein Aiolos antibody at 8µg/ml + HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 59, 52, 54, 30, 48, 50, 55, 42, 36, 33, 43, 31 kDa
 Observed band size: 59 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.