



glaB1 Antibody, Biotin conjugated

Product Code	CSB-PA406363HD01BOP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6KWM0
Immunogen	Recombinant Bacteroides vulgatus Alpha-1,3-galactosidase B protein (19-592AA)
Raised In	Rabbit
Species Reactivity	Bacteroides vulgatus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-1,3-galactosidase B (EC 3.2.1.n1) (EC 3.2.1.n2) (Exo-alpha-galactosidase B1) (EC 3.2.1.22), glaB1
Immunogen Species	Bacteroides vulgatus (strain ATCC 8482 / DSM 1447 / JCM 5826 / NBRC 14291 / NCTC 11154)
Research Area	Others
Target Names	glaB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.