



aroD Antibody, HRP conjugated

Product Code	CSB-PA326162HB01SWW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24670
Immunogen	Recombinant Salmonella typhi 3-dehydroquinase dehydratase protein (1-252AA)
Raised In	Rabbit
Species Reactivity	Salmonella typhi
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-dehydroquinase dehydratase (3-dehydroquinase) (EC 4.2.1.10) (Type I DHQase) (Type I dehydroquinase) (DHQ1), aroD
Immunogen Species	Salmonella typhi
Research Area	Others
Target Names	aroD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.