



ALTA1 Antibody

Product Code	CSB-PA303588LA01AZV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P79085
Immunogen	Recombinant Alternaria alternata Major allergen Alt a 1 protein (19-157AA)
Raised In	Rabbit
Species Reactivity	Alternaria alternata
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major allergen Alt a 1 (allergen Alt a 1), ALTA1
Immunogen Species	Alternaria alternata (Alternaria rot fungus) (Torula alternata)
Research Area	Others
Target Names	ALTA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.