



botA Antibody, Biotin conjugated

Product Code	CSB-PA320799LD01CLQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DPI0
Immunogen	Recombinant Clostridium botulinum Botulinum neurotoxin type A protein (1-436AA)
Raised In	Rabbit
Species Reactivity	Clostridium botulinum
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Botulinum neurotoxin type A (BoNT/A) (EC 3.4.24.69) (Bontoxilysin-A) (BOTOX) [Cleaved into: Botulinum neurotoxin A light chain; Botulinum neurotoxin A heavy chain], botA, atx bna
Immunogen Species	Clostridium botulinum
Research Area	Others
Target Names	botA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.