



# Outer capsid protein VP4 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA320838LB01ROH
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P11200
<b>Immunogen</b>	Recombinant Rotavirus A Outer capsid protein VP4 protein (247-479AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rotavirus A
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Outer capsid protein VP4 (Hemagglutinin) [Cleaved into: Outer capsid protein VP8*; Outer capsid protein VP5*]
<b>Immunogen Species</b>	Rotavirus A (strain RVA/Human/United Kingdom/ST3/1975/G4P2A[6]) (RV-A) (Rotavirus A (strain St. Thomas 3))
<b>Research Area</b>	Others
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.