



Outer capsid protein VP4 Antibody, Biotin conjugated

Product Code	CSB-PA320838LD01ROH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11200
Immunogen	Recombinant Rotavirus A Outer capsid protein VP4 protein (247-479AA)
Raised In	Rabbit
Species Reactivity	Rotavirus A
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer capsid protein VP4 (Hemagglutinin) [Cleaved into: Outer capsid protein VP8*; Outer capsid protein VP5*]
Immunogen Species	Rotavirus A (strain RVA/Human/United Kingdom/ST3/1975/G4P2A[6]) (RV-A) (Rotavirus A (strain St. Thomas 3))
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.