



DEFB126 Antibody

Product Code	CSB-PA006696EA01EXP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4H244
Immunogen	Recombinant Pongo pygmaeus Beta-defensin 126 protein (21-63AA)
Raised In	Rabbit
Species Reactivity	Pongo pygmaeus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-defensin 126 (Defensin, beta 126), DEFB126
Immunogen Species	Pongo pygmaeus (Bornean orangutan)
Research Area	Others
Target Names	DEFB126