



DEFB126 Antibody

Product Code	CSB-PA006696EA01EXP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4H244
Immunogen	Recombinant Pongo pygmaeus Beta-defensin 126 protein (21-63AA)
Raised In	Rabbit
Species Reactivity	Pongo pygmaeus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-defensin 126 (Defensin, beta 126), DEFB126
Immunogen Species	Pongo pygmaeus (Bornean orangutan)
Research Area	Others
Target Names	DEFB126
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.