



sacIM Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA523234LB01FNS |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O31073 |
| Immunogen | Recombinant Streptomyces achromogenes Modification methylase SacI protein (1-390AA) |
| Raised In | Rabbit |
| Species Reactivity | Streptomyces achromogenes |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Modification methylase SacI (M.SacI) (EC 2.1.1.37) (Cytosine-specific methyltransferase SacI), sacIM |
| Immunogen Species | Streptomyces achromogenes |
| Research Area | Others |
| Target Names | sacIM |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |