



sacIM Antibody, Biotin conjugated

Product Code	CSB-PA523234LD01FNS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O31073
Immunogen	Recombinant Streptomyces achromogenes Modification methylase SacI protein (1-390AA)
Raised In	Rabbit
Species Reactivity	Streptomyces achromogenes
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Modification methylase SacI (M.SacI) (EC 2.1.1.37) (Cytosine-specific methyltransferase SacI), sacIM
Immunogen Species	Streptomyces achromogenes
Research Area	Others
Target Names	sacIM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.