



LCY1 Antibody, HRP conjugated

Product Code	CSB-PA657129LB01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q38933
Immunogen	Recombinant Arabidopsis thaliana Lycopene beta cyclase, chloroplastic protein (81-501AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lycopene beta cyclase, chloroplastic (EC 5.5.1.19) (AtLCY) (Protein SUPPRESSOR OF ZEAXANTHIN-LESS 1), LCY1, LCYB LYC SZL1
Immunogen Species	Arabidopsis thaliana
Research Area	Others
Target Names	LCY1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.