



CD300LF Antibody, Biotin conjugated

Product Code	CSB-PA819470LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDQ1
Immunogen	Recombinant Human CMRF35-like molecule 1 protein (20-156AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CMRF35-like molecule 1 (CLM-1) (CD300 antigen-like family member F) (Immune receptor expressed on myeloid cells 1) (IREM-1) (Immunoglobulin superfamily member 13) (IgSF13) (NK inhibitory receptor) (CD antigen CD300f), CD300LF, CD300F CLM1 IGSF13 IREM1 NKIR
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD300LF