





HOXA1 Antibody

Product Code	CSB-PA010647DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49639
Immunogen	Recombinant Human Homeobox protein Hox-A1 protein (75-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Hox-A1 (Homeobox protein Hox-1F), HOXA1, HOX1F
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HOXA1


Image



Western blot
 All lanes: HOXA1 antibody at 2μg/ml
 Lane 1: NIH/3T3 whole cell lysate
 Lane 2: Hela whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 37, 15, 25 kDa
 Observed band size: 37 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA010647DSR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA010647DSR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.