





ICOSLG Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA010958ESR1HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O75144 |
| Immunogen | Recombinant Human ICOS ligand protein (19-256AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC; Recommended dilution: IHC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | ICOS ligand (B7 homolog 2) (B7-H2) (B7-like protein GI50) (B7-related protein 1) (B7RP-1) (CD antigen CD275), ICOSLG, B7H2 B7RP1 ICOSL KIAA0653 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | ICOSLG |

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA010958ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human epityphlon tissue using CSB-PA010958ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.