



POLR3A Antibody, HRP conjugated

Product Code	CSB-PA018343LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14802
Immunogen	Recombinant Human DNA-directed RNA polymerase III subunit RPC1 protein (392-632AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-directed RNA polymerase III subunit RPC1 (RNA polymerase III subunit C1) (EC 2.7.7.6) (DNA-directed RNA polymerase III largest subunit) (DNA-directed RNA polymerase III subunit A) (RNA polymerase III 155 kDa subunit) (RPC155) (RNA polymerase III subunit C160), POLR3A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLR3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.