



PLPPR2 Antibody, Biotin conjugated

Product Code	CSB-PA836219HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GM1
Immunogen	Recombinant Human Lipid phosphate phosphatase-related protein type 2 protein (1-156AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phospholipid phosphatase-related protein type 2 (EC 3.1.3.4) (Lipid phosphate phosphatase-related protein type 2) (Plasticity-related gene 4 protein) (PRG-4), PLPPR2, LPPR2 PRG4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLPPR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.