



Outer capsid glycoprotein VP7 Antibody, Biotin conjugated

Product Code	CSB-PA318053LD01ROH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10501
Immunogen	Recombinant Rotavirus A Outer capsid glycoprotein VP7 protein (51-326AA)
Raised In	Rabbit
Species Reactivity	Rotavirus A
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer capsid glycoprotein VP7
Immunogen Species	Rotavirus A
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.