



Reg3a Antibody

Product Code	CSB-PA019548LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35231
Immunogen	Recombinant Rat Regenerating islet-derived protein 3-alpha protein (26-174AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regenerating islet-derived protein 3-alpha (REG-3-alpha) (Islet of Langerhans regenerating protein 3) (REG 3) (Lithostathine 3) (Pancreatitis-associated protein 2) (RegIII) (Regenerating islet-derived protein III-alpha) (Reg III-alpha) [Cleaved into: Regenerating islet-derived protein 3-alpha 16.5 kDa form; Regenerating islet-derived protein 3-alpha 15 kDa form], Reg3a, Pap2 Reg3
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Reg3a
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.