



Prss29 Antibody

Product Code	CSB-PA859566LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99MS4
Immunogen	Recombinant Mouse Serine protease 29 protein (18-179AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease 29 (EC 3.4.21.-) (Implantation serine proteinase 2) (ISP-2) (Strypsin-2) (Strypsin-related protein) (Trypsin-like proteinase), Prss29, Isp2
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Prss29
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.