



Botulinum neurotoxin type E Antibody

Product Code	CSB-PA327013LA01DUK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30995
Immunogen	Recombinant Clostridium butyricum Botulinum neurotoxin type E protein (2-422AA)
Raised In	Rabbit
Species Reactivity	Clostridium butyricum
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Botulinum neurotoxin type E (BoNT/E) (EC 3.4.24.69) (Bontoxilysin-E) [Cleaved into: Botulinum neurotoxin E light chain; Botulinum neurotoxin E heavy chain]
Immunogen Species	Clostridium butyricum
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.