





PIAS3 Antibody

Product Code	CSB-PA896765DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6X2
Immunogen	Recombinant Human E3 SUMO-protein ligase PIAS3 protein (399-628AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 SUMO-protein ligase PIAS3 (EC 2.3.2.-) (E3 SUMO-protein transferase PIAS3) (Protein inhibitor of activated STAT protein 3), PIAS3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PIAS3

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA896765DSR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA896765DSR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.