





TNFRSF13B Antibody

Product Code	CSB-PA023971DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14836
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 13B protein (184-293AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 13B (Transmembrane activator and CAML interactor) (CD antigen CD267), TNFRSF13B, TACI
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TNFRSF13B

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA023971DSR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA023971DSR2HU at dilution of 1:100