





C1QB Antibody

Product Code	CSB-PA003638ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02746
Immunogen	Recombinant Human Complement C1q subcomponent subunit B protein (104-253AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement C1q subcomponent subunit B, C1QB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	C1QB

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA003638ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA003638ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.