



Recombinant Gallus gallus Syntaxin-6

Product Code	CSB-MP022902CH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZL19
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	EDPFFVVKGEVQKAVNTAQGLFQRWTELLQDPSIATREEIDWTTNELRNNLRS IEWDLEDLDETISIVEANPRKFNLDATELGIRKSFITSTRQVVREMKDQMSNSS MQALAERKNRQALLGESSQSWSGPDKYSRLDRELQLANSHFIEEQQAQQ QLIVEQQDEQLELVSGSIGVLKNMSQRIGGELEEQAVMLDDFSHELDSTHSRL DNVMKKLAKVSHMTSDR
Research Area	Others
Source	Mammalian cell
Target Names	STX6
Expression Region	4-231aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.