



Recombinant human Lipopolysaccharide-responsive and beige-like anchor protein

Product Code	CSB-EP013070HU(C)-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P50851
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	PQPHRHVLEISRQHEQPGQGIAPDAVNGQRRDSRSTVFRIFEFNWSQMHQRL LTDLLFSIETDIQMWRSHSTKTVMDVFNSSDNVIFVHNTIHLISQVMDNMVMAC GGILPLLSAATSATHELENIEPTQGLSIEASVTFLQRLISLVDVLIFASSLGFTEIE AEKSMSSGGILRQCLRLVCAVAVRNCLECCQHSQKTRGDKALKPMHSLIPLG KSAAKSPVDIVTGGISPV
Research Area	Signal Transduction
Source	E.coli
Target Names	LRBA
Expression Region	1267-1500aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.