



Recombinant human Lipopolysaccharide-responsive and beige-like anchor protein

| | |
|--------------------------|---|
| Product Code | CSB-YP013070HU(N) |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P50851 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | ASEDNRVPSPPTGDDGGGGGREETPTEGGALSLKPGLPIRGIRMKFAVLTGL VEVGEVSNRDIVETVFNLLVGGQFDLEMNFIQEGESINCMVDLLEKCDITCQA EVWSMFTAILKKSIRNLQVCTEVGLVEKVLGKIEKVDNMIADLLVDMLGVLASY NLTVRELKLFSSKLQGDKGRWPPHAGKLLSVLKHMPQKYGPDAFFNFP GKSA AAIALPPIAKWPYQNGFTFHTWLRMDPVNNINVDKDKP |
| Research Area | Signal Transduction |
| Source | Yeast |
| Target Names | LRBA |
| Expression Region | 1-251aa |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial of the full length of 2-2863aa |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |