



Recombinant human Testis-expressed sequence 12 protein

Product Code	CSB-MP023402HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BXU0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMANHLVKPDNRNCKRPRELESPVPDSPQLSSLGKSDSSFSEISGLFYKDEAL EKDLNDVSKEINLMLSTYAKLLSERA AVDASYIDEIDELFKEANAIENFLIQKREF LRQRFTVIANTLHR
Research Area	Cell Biology
Source	Mammalian cell
Target Names	TEX12
Expression Region	1-123aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Target Details	This gene is similar to a mouse gene that is expressed in the testis.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.