



Recombinant human Protein FAM110B

| | |
|--------------------------|---|
| Product Code | CSB-YP008005HU |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8TC76 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MPTETLQTGSMVKPVSPAGTFTSAVPLRILNKGPDYFRRQAEPNPKRLSAVER LEADKAKYVKSQEVINAKQEPVKPAVLAKPPVCPAAKRALGSPTLKVFGNHAK TESGVQRENKLEILKNIINSSEGSSSGSGHGHSSRNWPPHRSEATDLHRHSF AESLKVYPTQGRRSPQEGGSHVGRRLLEQSAESFLHVSHSSSDIRKVTSVKPL KAIPCSSSAPPLPPKPKIAAIIASMSPEADPVEPACGVSRPRLQRSKSDLSDR YFRVDADVERFFNYCGLDPEELENLGMENFARANSIIISLNFRRSASMISSDCE QSQDSNSDLRNDSDSANDRVYPYGISAIIERNARIKWLYSIKQARESQKVSHV |
| Research Area | Cancer |
| Source | Yeast |
| Target Names | FAM110B |
| Expression Region | 1-370aa |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |