



Recombinant human Protein FAM110B

Product Code	CSB-YP008005HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TC76
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPTETLQTGSMVKPVSPAGTFTSAVPLRILNKGPDYFRRQAEPNPKRLSAVER LEADKAKYVKSQEVINAKQEPVKPAVLAKPPVCPAAKRALGSPTLKVFGNHAK TESGVQRENKLEILKNIINSSEGSSSGSGHGHSSRNWPPHRSEATDLHRHSF AESLKVYPTQGRRSPQEGGSHVGRRLLEQSAESFLHVSHSSSDIRKVTSVKPL KAIPCSSSAPPLPPKPKIAAIASMKSP EADPVEPACGVSRPSLQRSKSDLSDR YFRVDADVERFFNYCGLDPEELENLGMENFARANSDIISLNFRRSASMISSDCE QSQDSNSDLRNDSDANDRVYPYGISA IERNARIKWLYSIKQARES QKVSHV
Research Area	Cancer
Source	Yeast
Target Names	FAM110B
Expression Region	1-370aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.