




Recombinant Rabbit Vascular endothelial growth factor A (VEGFA) (X303W)

Product Code	CSB-YP025833RB
Abbreviation	Recombinant Rabbit VEGFA protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	XP_002714743.1
Product Type	Recombinant Proteins
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	YLWLLHAPLLHVSLDSLVAAFSVVFEVLPSSLDIPGSKKEEQASTQLGGVEEEH EASGVVTRQLAFVVDAITNPTLEKETKTITRIHKAPKFPSHFGFWKRAEEERGG KSSGEKSRKTDGVGERAQAGEQREGPGQRAAALTDRTDTAPSPSAHLLPG RRPTVDAAASGGQQPEPAPEGGVGARGIALKLFVQLLGCSRSGGVVRA GGAEPSCAARSVSSGREPPPPPEEEGEEGEEKEEERGPWRGLGAGEPGS WTGEAAVCADSAPAARAPQALARASAPGARGARGGAEESGPSRSPSRRGSA SRAGPGRASETMNFLLSVHWSLALLLYLHHAKWSQAAPMAEEGDNKPHEV VKFMEVYRRSYCQPIETLVDIFQEYPDEIEYIFKPSCVPLVRCGGCCNDESLEC VPTEEFNVTMIMRIKPHQGGQHIGEMSFLQHNKCECRCDKPRR
Research Area	Others
Source	Yeast
Target Names	VEGFA
Expression Region	51-511aa (X303W)
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	51.7kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.