



# Recombinant human Transcription initiation factor IIB

<b>Product Code</b>	CSB-BP010001HU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q00403
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASTSRLDALPRVTCPNHPDAILVEDYRAGDMICPECGLVVGDRVIDVGSEWR TFSNDKATKDP SRV GDSQNPLLS DGD LSTMIGKGTGAASFDEF GNSKYQNR TMSSSDRAMMNAFKEITTMADRINLPRNIVDR TNNL FKQVYE QKSLKGRAN IASACLYIACRQEGVPRTFKEICAVSRISKKEIGRCFKLILKALET SVDLIT TGDFM SRFC SNLCLPKQVQMAATHIARKAVELDLPGRSPISVAAAAIYMASQ ASAEKR TQKEIGDIAGVADV TIRQSYRLIYPRAPDLFPTDFKFDTPVDKLPQL
<b>Research Area</b>	Immunology
<b>Source</b>	Baculovirus
<b>Target Names</b>	GTF2B
<b>Expression Region</b>	1-316aa
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length
<b>Target Details</b>	This gene encodes the general transcription factor IIB, one of the ubiquitous factors required for transcription initiation by RNA polymerase II. The protein localizes to the nucleus where it forms a complex (the DAB complex) with transcription factors IID and IIA. Transcription factor IIB serves as a bridge between IID, the factor which initially recognizes the promoter sequence, and RNA polymerase II.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.