



Recombinant human F-box/LRR-repeat protein 18

Product Code	CSB-EP008469HU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96ME1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MASSGEDISNDDDDMHPAAAGMADGVHLLGFSDEILLHILSHVPSTDLILNVRR TCRKLAALCLDKSLIHTVLLQKDYQASEDKVRQLVKEIGREIQQLSMAGCYWLP GSTVEHVARCRSLVKVNLGCHLTSRLSKMLSALQHLSLAIDVSPGFDASQ LSSECKATLSRVRELKQTLFTPSYGVVPCCTSLEKLLLYFEILDRTREGAILSGQ LMVGQSNVPHYQNLRVFYARLAPGYINQEVVRLYLAVLSDRTPQNLHAFLISV PGSFAESGATKNLLDSMARNVVDALQLPKSWLNGSSLLQHMKFNPFYFSF SRCTLSGGHLIQQVINGGKDLRSLASLNLSGCVHCLSPDSLLCR
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	FBXL18
Expression Region	1-365aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Isoform3
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.