



Recombinant Human DNA-binding protein inhibitor ID-2

Product Code	CSB-YP010967HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q02363
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKAFSPVRSVRKNSLSDHSLGISRSKTPVDDPMSLLYNMND CYSK LKELVPSI PQNKKVSKMEILQHVIDYILD LQIALDSHPTIVSLHHQRPGQNQASRTPLTTLNT DISILSLQASEFPSELMSNDSKALCG
Research Area	Epigenetics and Nuclear Signaling
Source	Yeast
Target Names	ID2
Expression Region	1-134aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Target Details	This protein belongs to the inhibitor of DNA binding (ID) family, members of which are transcriptional regulators that contain a helix-loop-helix (HLH) domain but not a basic domain. Members of the ID family inhibit the functions of basic helix-loop-helix transcription factors in a dominant-negative manner by suppressing their heterodimerization partners through the HLH domains. This protein may play a role in negatively regulating cell differentiation. A pseudogene has been identified for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.