



Recombinant Bacillus subtilis (strain 168) Inositol 2-dehydrogenase/D-chiro-inos

Product Code	CSB-YP329675BRG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P26935
Product Type	Recombinant Protein
Immunogen Species	Bacillus subtilis (strain 168)
Purity	≥85% (SDS-PAGE)
Sequence	MSLRIGVIGTGAIGKEHINRITNKLSGAEIVAVTDVNQEEAAQKVVEQYQLNATVY PNDDSLADENVDAVLVTSWGPAHESSVLKAIKAQKYVFCEKPLATTAEGCMR IVEEEIKVGRKRLVQVGFMRRYDSGYVQLKEALDNHVICEPLMIHCAHRNPTVG DNYTTDMAVVDLTVHEIDVLHVLNDDYESVQVIYPKSKNALPHLKDPQIVVI ETKGGIVINAEIYVNCKYGYDIQCEIVGEDGIIKLPEPSSISLRKEGRFSTDILMD WQRRFVAAYDVEIQDFIDSIQKKGEVSGPTAWDGYIAAVTTDACVKAQESGQK EKVELKEKPEFYQSFTTVQN
Research Area	Others
Source	Yeast
Target Names	iolG
Expression Region	1-344aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.