



Recombinant Epstein-Barr virus Deoxyuridine 5'-triphosphate nucleotidohydrolase (BLLF3)

Product Code	CSB-MP319032EFC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	K9US42
Product Type	Recombinant Protein
Purity	≥85% (SDS-PAGE)
Sequence	MEACPHIRYAFQNDKLLLQQASVGRLLTLVNKTTILLRPMKTTTVDLGLYARPPE GHGLMLWGSTSRPVTSHVGIIDPGYTGELRLILQNQRRYNSTLRPSELKIHAA FRYATPQMEEDKGPINHPQYPGDVGLDVSLPKDLALFPHQTVSVTLTVPPPSI PHHRPTIFGRSGLAMQGILVKPCRWRGGVDVSLTNFSDQTVFLNKYRRFCQ LVYLHKHHLTSFYSPHSDAGVLGPRSLFRWASCAFEVPSLAMGDSGLSEAL EGRQGRGFGSSGQ
Research Area	Others
Source	Mammalian cell
Expression Region	1-278aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.