



Recombinant Human Heterogeneous nuclear ribonucleoprotein A1 (HNRNPA1), partial

Product Code	CSB-YP010600HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09651
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SKSESPKEPEQLRKLFIGGLSFETTDESLRSHFEQWGLTDCVVMRDPNTRKRS RGFGFVTYATVEEVDAAMNARPHKVDGRVVEPKRAVSREDSQRPGAHLTVK KIFVGGIKEDTEEHHLRDYFEQYGKIEVIEIMDRGSGKKRGFAFVTFDDHDSV DKIVIQKYHTVNGHNCEVRKALSKQEMASASSSQRGRSGSGNFGGGRGGGF GGNDNFGRGGNFSGRGGFGGSRGGGGYGGSGDGYNGFGNDGGYGGGGP GYSGGSRGYGSQQGYGNQGSYGGSGSYDSYNNGGGGGFGGGSGSNF GGGGSYNDFGNYNNQSSNFGPMKGGNFGRSSGPYGGGGQYFAKPRNQ
Source	Yeast
Target Names	HNRNPA1
Protein Names	Recommended name: Heterogeneous nuclear ribonucleoprotein A1 Short name= hnRNP A1 Alternative name(s): Helix-destabilizing protein Single-strand RNA-binding protein hnRNP core protein A1
Expression Region	2-354aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.