



# Recombinant Mouse Ccn family member 2 (Ccn2), partial

<b>Product Code</b>	CSB-MP006147MO
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P29268
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	QDCSAQCQCAAEAAPHCPAGVSLVLDGCGCCRVCAKQLGELCTERDPCDPH KGLFCDFGSPANRIGVCTAKDGAPCVFGGSVYRSGESFQSSCKYQCTCLDG AVGCVPLCSMDVRLPSPDCPFPRRVKLPGKCCEEWVCDEPKDRTAVGPALA AYRLEDTFGPDPTMMRANCLVQTTEWSACSKTCGMGISTRVTNDNTFCRLEK QSRLCMVRPCEADLEENIKKGGKCI RTPKIAKPVKFE LSGCTSVKTYRAKFCGV CTDGRCTPHRTTTL PVEFKCPDGEIMKKNMMFIKTCACHYNC PGDNDIFESL YYRKMY
<b>Research Area</b>	Others
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Ccn2
<b>Expression Region</b>	26-344aa
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Target Details</b>	This protein is a mitogen that is secreted by vascular endothelial cells. The encoded protein plays a role in chondrocyte proliferation and differentiation, cell adhesion in many cell types, and is related to platelet-derived growth factor. Certain polymorphisms in this gene have been linked with a higher incidence of systemic sclerosis.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.