



Recombinant Human Homeobox protein OTX2 (OTX2)

Product Code	CSB-MP017299HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32243
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMSYLKQPPYAVNGLSLTTS GMDLLHPSVGYGPGWASCPAATPRKQRRERT TFTRAQLDVLEALFAKTRYPDIFMREEVALKINLPESRVQVWFKNRRRAKCRQQ QQQQQNGGQNKVRPAKKKTSPAREVSSSESGTSGQFTPPSSTSVPTIASSSAP VSIWSPASISPLSDPLSTSSSCMQRSYPMTYTQASGYSQGYAGSTSYFGGMD CGSYLTPMHHQLPGPGATLSPMG TNAV TSHLNQSPASLSTQGYGASSLGFNS TTDCLDYKDQTASWKLNFNADCLDYKDQTSSWK FQVL
Research Area	Epigenetics and Nuclear Signaling
Source	Mammalian cell
Target Names	OTX2
Expression Region	1-297aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Isoform2
Target Details	This gene encodes a member of the bicoid sub-family of homeodomain-containing transcription factors. The encoded protein acts as a transcription factor and may play a role in brain and sensory organ development. A similar protein in mice is required for proper forebrain development. Two transcript variants encoding distinct isoforms have been identified for this gene. Other alternative splice variants may exist, but their full length sequences have not been determined.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.