



Recombinant Cricetulus griseus Cathepsin D (Ctsd), partial

Product Code	CSB-BP006187DXU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	G3I4W7
Product Type	Recombinant Protein
Purity	≥85% (SDS-PAGE)
Sequence	GPVSELLKNYLDAQYYGEIGIGTTPQCFTVVFDTGSSNLWVPSIHCKLLDIACW IHHKYNSGKSSTFVKNGTSFDIHYGSGLSGYLSQDTVSVPCCKSEQPGGLKVE KQIFGEAIKQPGITFIAAKFDGILGMGYPSISVNNVVPVFDNLMQQKLVEKNIFS FLNRDPTGQPGGELMLGGIDSKYYEGELSYLNVTRKAYWQVHMDQLDVANG LTLCKGGCEAIVDTGTSLLVGPVDEVKELQKAIGAVPLIQGEYMIPCEKVSSLPS VTLKLGKDYELSPSKYVLKVSQGGKTICLSGFMGMDIPPPSGPLWILGDVFIG TYYTVFDRDNNRVGFAKAATL
Source	Baculovirus
Expression Region	65-408aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.