



# Recombinant Cricetulus griseus Cathepsin D (Ctsd), partial

<b>Product Code</b>	CSB-EP006187DXU-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	G3I4W7
<b>Product Type</b>	Recombinant Protein
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GPVSELLKNYLDAQYYGEIGIGTTPQCFTVVFDTGSSNLWVPSIHCKLLDIACW IHHKYNSGKSSTFVKNGTSFDIHYGSGLSGYSQDTVSVPCCKSEQPGGLKVE KQIFGEAIKQPGITFIAAKFDGILGMGYPSISVNNVVPVFDNLMQQKLVEKNIFS FLNRDPTGQPGGELMLGGIDSKYYEGELSYLVNTRKAYWQVHMDQLDVANG LTLCKGGCEAIVDTGTSLLVGPVDEVKELQKAIGAVPLIQGEYMIPCEKVSSLPS VTLKLGKDYELSPSKYVLKVSQGGKTICLSGFMGMDIPPPSGPLWILGDVFIG TYYTVFDRDNNRVGFAKAATL
<b>Source</b>	E.coli
<b>Expression Region</b>	65-408aa
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.