



IL18BP Monoclonal Antibody

Product Code	CSB-MA222469
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human IL18BP protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of IL18BP and does not cross-react with related proteins.
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Relevance	Isoform A binds to IL-18 and inhibits its activity.Functions as an inhibitor of the early TH1 cytokine response.
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG2b
Clonality	Monoclonal
Alias	I18BP_HUMAN; IL-18BP; IL18 BP; IL18 BPα; IL18BP; IL18BPα; Interleukin 18 binding protein; Interleukin-18-binding protein; MC51L 53L 54L homolog gene product; Tadekinig alfa; Tadekinig-alfa.
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	12C2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.