



PPP1R18 Monoclonal Antibody

Product Code	CSB-MA940287
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human PPP1R18 protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of PPP1R18 and does not cross-react with related proteins.
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Relevance	May target protein phosphatase 1 to F-actin cytoskeleton.
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Clonality	Monoclonal
Alias	HKMT1098; KIAA1949; Protein phosphatase 1 F actin cytoskeleton targeting subunit; PPP1R18;
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)