



SMYD4 Monoclonal Antibody

Product Code	CSB-MA278216
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human SMYD4 protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of SMYD4 and does not cross-react with related proteins.
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Clonality	Monoclonal
Alias	KIAA1936; SET and MYND domain containing 4; SET and MYND domain-containing protein 4; SMYD4; SMYD4_HUMAN; ZMYND21.
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)