



SMYD4 Monoclonal Antibody

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-MA278216 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Immunogen | Purified recombinant Human SMYD4 protein fragments expressed in E.coli. |
| Raised In | Mouse |
| Species Reactivity | Human |
| Specificity | This antibody detects endogenous levels of SMYD4 and does not cross-react with related proteins. |
| Tested Applications | ELISA,WB;Recommended dilution:WB:1:500-1:5000 |
| Form | Ascites |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Alias | KIAA1936; SET and MYND domain containing 4; SET and MYND domain-containing protein 4; SMYD4; SMYD4_HUMAN; ZMYND21. |
| Product Type | Monoclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |