



GATA6 Monoclonal Antibody

Product Code	CSB-MA869500
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human GATA6 protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Relevance	Thought to be important for regulating terminal differentiation and/or proliferation.
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	unknown
Clonality	Monoclonal
Alias	ASD9, AVSD5, PACHD
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	3H3-H6-H11