



GATA6 Monoclonal Antibody

Product Code	CSB-MA869500
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human GATA6 protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Relevance	Thought to be important for regulating terminal differentiation and/or proliferation.
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	unknown
Clonality	Monoclonal
Alias	ASD9, AVSD5, PACHD
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	13F4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.