



# GATA6 Monoclonal Antibody

<b>Product Code</b>	CSB-MA916172
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Purified recombinant Human GATA6 protein fragments expressed in E.coli.
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	E WB Recommended dilution: Predicted MW: 60kd, Western blotting: 1:1000
<b>Relevance</b>	Thought to be important for regulating terminal differentiation and/or proliferation.
<b>Form</b>	Ascites
<b>Isotype</b>	unknown
<b>Clonality</b>	Monoclonal
<b>Alias</b>	ASD9, AVSD5, PACHD
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Clone No.</b>	3H3-H6-H11