



SMYD1 Monoclonal Antibody

Product Code	CSB-MA898680
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human SMYD1 protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Relevance	Acts as a transcriptional repressor. Essential for cardiomyocyte differentiation and cardiac morphogenesis By similarity.
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Clonality	Monoclonal
Alias	BOP, KMT3D, ZMYND18, ZMYND22
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)