



TXNRD2 Monoclonal Antibody

Product Code	CSB-MA576515
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human TXNRD2 protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Relevance	Maintains thioredoxin in a reduced state. Implicated in the defenses against oxidative stress. May play a role in redox-regulated cell signaling. Maintains thioredoxin in a reduced state. Implicated in the defenses against oxidative stress.
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Clonality	Monoclonal
Alias	TR, TR3, SELZ, TRXR2, TR-BETA
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	18F5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.